



# LookinBody Web

## InBody LAN Connection Guide



# Introduction

---

Compatible InBody units can connect via a LAN/Ethernet connection. A physical wired connection provides a more stable and secure connection to your facility's network. Follow this guide to do so.

## **Needed:**

- Cat5/Cat5e/Cat6 Cable
- Administrator access to the InBody

**Estimated Time to Complete: 10 Minutes**

# Connecting via Ethernet Cable/LAN


---



# Setting Up the LAN Connection

---

Connect the ethernet cable to the InBody unit




# Setting Up the LAN Connection

---



Tap the Administrator Menu button.


# Setting Up the LAN Connection



Select 16. Internet Options


# Setting Up the LAN Connection

---



Tap the O to enable connection.

# Setting Up the LAN Connection




Tap the check mark box to the left of LAN (Wired)

Tap Next.




# Setting Up the LAN Connection




Tap Next.

Please wait, while the InBody finalizes the setup




# Setting Up the LAN Connection



Confirm that that a valid IP Address is displayed.

Tap **Exit** to leave this screen.

# Device Login



After connecting your unit via LAN,  
select **04. Cloud Service**

# Device Login

**04. Cloud Service** Help Exit

If you enable the 'Cloud Service', members will have the option to check their InBody results from the registered mobile device.

Do you want to enable the Cloud Service?

Lookin'Body Website Account Login

You can manage the InBody results from the web. To create an account, please contact Customer Service.

The contact information for Customer Service is shown below.

**InBody**  
TEL : 1-323-932-6503

Save

1. Tap O to enable the service
2. Tap the Check Mark Box to open the Login window

# Device Login

**Lookin'Body Website Account Login** Exit

**ID**

**Password**

If you forgot your password, please contact Customer Service.

The contact information for Customer Service is shown below.  
Biospace, Inc.  
TEL : 1-323-932-6503

**Login**

Tap on each field to update the login ID and password with the information below

Your Device ID and Device Password can be found in the account registration email you received when you first signed up for your account.

# Device Login

---



Lookin'Body Website Account Login Exit

ID  
InBody

Password  
\*\*\*\*\*

If you forgot your password, please contact Customer Service.

The contact information for Customer Service is shown below.  
Biospace, Inc.  
TEL : 1-323-932-6503

Login

Tap Login.

If the unit displays an error message, contact our support team at:

323-932-6503 extension 2

*or*

[LBWeb@InBody.com](mailto:LBWeb@InBody.com)

# Device Login



**04. Cloud Service** Help Exit

If you enable the 'Cloud Service', members will have the option to check their InBody results from the registered mobile device.

Do you want to enable the Cloud Service? ✕ ○

☒ **Lookin'Body Website Account Login**  
Already logged in.

**Save**

If successful, the Check Mark Box will be highlighted in red.

If permitted, tap **Save**. If not tap **Exit** on the top right.

# Verify Unit Connection Status




Successfully connected units will display a cloud icon with “LB”



# Verify Unit Connection Status

Depending on the version of the unit, the main screen will display 2-3 different types of status icons:



— If the top right corner does not display a cloud symbol, the InBody unit is disconnected from the network and internet. Please review pages 6-18.



— If the top right corner displays a cloud symbol, the unit is connected to the internet and ready to go. If the unit was restarted or the personal profile information does not appear automatically, log the unit back in as shown in slides 12-17.



— In new versions of the InBody unit, if a cloud symbol with the letters “LB” appears, the unit is connected and logged in. No further action is required.

# Additional Information

---



## Overview

- ❗ Compatible InBody units can connect via the internet using either a Wi-Fi connection or LAN connection.
- ❗ A LAN connection provides a more stable and secure connection to the internet; however, it may be limited according to the length of the cable and available ports on the router.
- ❗ If your computer network is managed by an IT team, it is recommended to:
  - Speak to the IT team before attempting connection as some networks may require additional security precautions.

**For any issues or questions, please feel free to contact the support team for assistance.**

---

**LookinBody Web  
Support Team**



[LBWeb@InBody.com](mailto:LBWeb@InBody.com)



323.932.6503 extension 2



[www.InBody.com](http://www.InBody.com)

